

Indication Matrix – Hand

Mode of action	Stabilisation							Immobilisation	
	slight			moderate		high			
	Diagonal Thumb Support	Manu Rheuma Flex	Manu Sensa	Rhizo Forsa	Manu Comfort	Thumboform Long	Manu Forsa Volar	Manu Immobil	Manu Immobil Long
Article number	4074	4024	50P13	50A4	4055	4085	50P12	50P10	50P11
Indication									
Thumb									
Irritation of the thumb	•			•		•			
Osteoarthritis of the thumb saddle joint	•			•		•			
Ulnar collateral ligament injury (skier's thumb)	•			•		•			
De Quervain's tenosynovitis	•					•			
Excessive strain	•								
Wrist joint									
Rheumatoid arthritis		•	•						
Osteoarthritis of the wrist		•	•		•	•	•	•	•
Excessive strain		•					•		
Irritation of the wrist		•	•		•		•	•	•
Tendosynovitis		•	•		•		•		
Feeling of instability			•						
Carpal tunnel syndrome			•		•		•	•	•
Wrist sprain			•		•		•	•	•
Wrist dislocation					•		•	•	•
Posttraumatic and postoperative treatment		•	•		•		•	•	•
Paresis/paralysis					•		•	•	•
Hemiplegia					•				•
Postoperative positioning, e.g. after wrist fractures								•	•

Manu ComforT

Wrist Support

The 4055 Manu ComforT prevents improper wrist movement and provides reliable stabilisation. Flexion and extension of the wrist are limited by an anatomically pre-shaped volar splint, dorsal coil springs and a central stabilising strap.




Benefits

- Skin-friendly, dirt resistant material ensures a high degree of wearing comfort
- Individually, continuously adjustable straps allow for a proper fit
- Strap guide loops support the quick and easy handling

Prescription note

Wrist support with adjustable compression of soft tissue and a volar splint to stabilise the wrist



PCM technology

The soft, temperature-regulating material PCM (Phase Change Material) ensures good wearing comfort as it equalises temperature fluctuations. Micro-pellets integrated in the material melt or solidify to achieve a cooling or warming effect. This unique function is maintained even after long periods of use and repeated laundering.

Indications

- Acute and chronic irritation
- Functional treatment of tendosynovitis
- Carpal tunnel syndrome
- Wrist sprain
- Wrist dislocation
- Posttraumatic and postoperative treatment
- Osteoarthritis of the wrist
- Paresis/paralysis
- Hemiplegia

Quality for life